


## SHIGEMI - NMR Pricelist Version 01.01.2024



## Index

Contents	Page
<b>SHIGEMI NMR TUBES</b>	
SHIGEMI TUBE SET for Acetone D6	3
SHIGEMI TUBE SET for D2O	4
SHIGEMI TUBE SET for CDCl <sub>3</sub>	5
SHIGEMI TUBE SET DMSO-D6	6
SHIGEMI TUBE SET Methanol-D4	7
SHIGEMI Precision Large Volume NMR Tubes	8
SHIGEMI Precision Large Volume NMR Tubes	8
SHIGEMI Precision Coaxial System	9
SHIGEMI 5mm Micro Bottom Tubes	10
SHIGEMI Quartz NMR Tubes	11
SHIGEMI Quartz ESR (EPR) Tubes	11
SHIGEMI Slot Tube	11
DURAN® 5mm Ultra Precision Thin Wall NMR Sample Tubes	12
Pyrex® 4mm Ultra Precision Thin Wall NMR Sample Tubes	12
Pyrex® 3mm Ultra Precision Thin Wall NMR Sample Tubes	12
Our Partners	13

## SHIGEMI TUBES

### SHIGEMI TUBE SET for Acetone D6

Article Number	Description	Quantity	Price
AMS-0025	2,5mm Shigemi TUBE SET for Acetone-D6 for BRUKER Outer OD 2.5mm, Insert OD 1.8mm, Insert length 115mm, Outer length	1	190,00 €
AMS-003	3mm Shigemi TUBE SET for Acetone-D6 for every System Outer OD 3.0mm, Insert OD 2.4mm, Insert length 190mm, Outer length	1	160,00 €
AMS-005B	5mm Shigemi TUBE SET for Acetone-D6 for BRUKER Outer OD 5.0mm, Insert OD 4.1mm, Insert length 190mm, Outer length 180mm, 8mm Bottom length	1	165,00 €
AMS-005J	5mm Shigemi TUBE SET for Acetone-D6 for JOEL Outer OD 5.0mm, Insert OD 4.1mm, Insert length 190mm, Outer length	1	160,00 €
AMS-005V	5mm Shigemi TUBE SET for Acetone-D6 for VARIAN Outer OD 5.0mm, Insert OD 4.1mm, Insert length 190mm, Outer length	1	160,00 €
AMS-008BJ	8mm Shigemi TUBE SET for Acetone-D6 for BRUKER and JOEL Outer OD 8.0mm, Insert OD 6.9mm, Insert length 200mm, Outer length 190mm, 10mm Bottom length	1	410,00 €
AMS-008V	8mm Shigemi TUBE SET for Acetone-D6 for VARIAN Outer OD 8.0mm, Insert OD 6.9mm, Insert length 200mm, Outer length 190mm, 15mm Bottom length	1	410,00 €
AMS-010BJ	10mm Shigemi TUBE SET for Acetone-D6 for BRUKER and JOEL Outer OD 10.0mm, Insert OD 8.9mm, Insert length 200mm, Outer length 190mm, 10mm	1	460,00 €
AMS-010V	10mm Shigemi TUBE SET for Acetone-D6 for VARIAN Outer OD 10.0mm, Insert OD 8.9mm, Insert length 200mm, Outer length 190mm, 15mm Bottom length	1	460,00 €

### SHIGEMI TUBE SET for D2O

Article Number	Description	Quantity	Price
----------------	-------------	----------	-------

BMS-0025	2,5mm Shigemi TUBE SET for D2O for BRUKER Outer OD 2.5mm, Insert OD 1.8mm, Insert length 115mm, Outer length 120mm, 6mm Bottom length	1	195,00 €
BMS-003	3mm Shigemi TUBE SET for D2O for every System Outer OD 3mm, Insert OD 2.4mm, Insert length 190mm, Outer length 180mm, 10mm Bottom	1	155,00 €
BMS-004B	4mm Shigemi TUBE SET for D2O for BRUKER Outer OD 4mm, Insert OD 3.1mm, Insert length 190mm, Outer length 180mm, 8mm Bottom	1	155,00 €
BMS-004J	4mm Shigemi TUBE SET for D2O for JEOL Outer OD 4mm, Insert OD 3.1mm, Insert length 190mm, Outer length 180mm, 12mm Bottom	1	155,00 €
BMS-004V	4mm Shigemi TUBE SET for D2O for VARIAN Outer OD 4mm, Insert OD 3.1mm, Insert length 190mm, Outer length 180mm, 15mm Bottom	1	155,00 €
BMS-005B	5mm Shigemi TUBE SET for D2O for BRUKER Outer OD 5mm, Insert OD 4.1mm, Insert length 190mm, Outer length 180mm, 8mm Bottom	1	120,00 €
BMS-005J	5mm Shigemi TUBE SET for D2O for JEOL Outer OD 5mm, Insert OD 4.1mm, Insert length 190mm, Outer length 180mm, 12mm Bottom	1	140,00 €
BMS-005V	5mm Shigemi TUBE SET for D2O for VARIAN Outer OD 5mm, Insert OD 4.1mm, Insert length 190mm, Outer length 180mm, 15mm Bottom length	1	140,00 €
BMS-008BJ	8mm Shigemi TUBE SET for D2O for BRUKER and JOEL Outer OD 8mm, Insert OD 6.9mm, Insert length 200mm, Outer length 190mm, 10mm Bottom	1	420,00 €
BMS-008V	8mm Shigemi TUBE SET for D2O for VARIAN Outer OD 8mm, Insert OD 6.9mm, Insert length 200mm, Outer length 190mm, 15mm Bottom	1	420,00 €
BMS-010BJ	10mm Shigemi TUBE SET for D2O for BRUKER and JOEL Outer OD 10mm, Insert OD 8.9mm, Insert length 200mm, Outer length 190mm, 10mm Bottom length	1	490,00 €
BMS-010V	10mm Shigemi TUBE SET for D2O for VARIAN Outer OD 10mm, Insert OD 8.9mm, Insert length 200mm, Outer length 190mm, 15mm Bottom length	1	545,00 €
BMS-005TB	5mm Shigemi Thin Wall TUBE SET ID 4.52mm for D2O for BRUKER Outer OD 5mm, Insert OD 4.1mm, Insert length 190mm, Outer length 180mm, 8mm Bottom	1	170,00 €
BMS-005TJ	5mm Shigemi Thin Wall ID 4.52mm TUBE SET for D2O for JOEL Outer OD 5mm, Insert OD 4.1mm, Insert length 190mm, Outer length 180mm, 12mm Bottom	1	170,00 €
BMS-005TV	5mm Shigemi Thin Wall ID 4.52mm TUBE SET for D2O for VARIAN Outer OD 5mm, Insert OD 4.1mm, Insert length 190mm, Outer length 180mm, 15mm Bottom length	1	170,00 €
BMS-008TB	8mm Shigemi Thin Wall TUBE SET ID 7.5mm for D2O for BRUKER Outer OD 8mm, Insert OD 6.9mm, Insert length 200mm, Outer length 190mm, 10mm Bottom length	1	440,00 €
BMS-008TV	8mm Shigemi Thin Wall TUBE SET ID 7.5mm for D2O for VARIAN Outer OD 8mm, Insert OD 6.9mm, Insert length 200mm, Outer length 190mm, 15mm Bottom length	1	440,00 €

### SHIGEMI TUBE SET for CDCI3

Article Number	Description	Quantity	Price
CMS-0025	2,5mm Shigemi TUBE SET for CDCI3 for BRUKER Outer OD 2.5mm, Insert OD 1.8mm, Insert length 115mm, Outer length 120mm, 6mm Bottom	1	185,00 €
CMS-003	3mm Shigemi TUBE SET for CDCI3 for every System Outer OD 3.0mm, Insert OD 2.4mm, Insert length 190mm, Outer length 180mm, 10mm Bottom	1	155,00 €
CMS-005B	5mm Shigemi TUBE SET for CDCI3 for BRUKER Outer OD 5.0mm, Insert OD 4.1mm, Insert length 190mm, Outer length 180mm, 8mm Bottom length	1	120,00 €
CMS-005J	5mm Shigemi TUBE SET for CDCI3 for JOEL Outer OD 5.0mm, Insert OD 4.1mm, Insert length 190mm, Outer length 180mm, 12mm Bottom	1	140,00 €
CMS-005V	5mm Shigemi TUBE SET for CDCI3 for VARIAN Outer OD 5.0mm, Insert OD 4.1mm, Insert length 190mm, Outer length 180mm, 15mm Bottom	1	140,00 €
CMS-008BJ	8mm Shigemi TUBE SET for CDCI3 for BRUKER and JOEL Outer OD 8.0mm, Insert OD 6.9mm, Insert length 200mm, Outer length 190mm, 10mm Bottom length	1	420,00 €
CMS-008V	8mm Shigemi TUBE SET for CDCI3 for VARIAN Outer OD 8.0mm, Insert OD 6.9mm, Insert length 200mm, Outer length 190mm, 15mm Bottom	1	420,00 €
CMS-010BJ	10mm Shigemi TUBE SET for CDCI3 for BRUKER and JOEL Outer OD 10.0mm, Insert OD 8.9mm, Insert length 200mm, Outer length 190mm, 10mm	1	490,00 €
CMS-010V	10mm Shigemi TUBE SET for CDCI3 for VARIAN Outer OD 10.0mm, Insert OD 8.9mm, Insert length 200mm, Outer length 190mm, 15mm Bottom length	1	540,00 €
CMS-005TB	5mm Shigemi Thin Wall TUBE SET ID 4.52mm for CDCI3 for BRUKER Outer OD 5mm, Insert OD 4.1mm, Insert length 190mm, Outer length 180mm, 8mm Bottom	1	170,00 €
CMS-005TJ	5mm Shigemi Thin Wall ID 4.52mm TUBE SET for CDCI3 for JOEL Outer OD 5mm, Insert OD 4.1mm, Insert length 190mm, Outer length 180mm, 12mm Bottom	1	170,00 €
CMS-005TV	5mm Shigemi Thin Wall ID 4.52mm TUBE SET for CDCI3 for VARIAN Outer OD 5mm, Insert OD 4.1mm, Insert length 190mm, Outer length 180mm, 15mm Bottom length	1	170,00 €
CMS-008TB	8mm Shigemi Thin Wall TUBE SET ID 7.5mm for D2O for BRUKER Outer OD 8mm, Insert OD 6.9mm, Insert length 200mm, Outer length 190mm, 10mm Bottom length	1	440,00 €
CMS-008TV	8mm Shigemi Thin Wall TUBE SET ID 7.5mm for D2O for VARIAN Outer OD 8mm, Insert OD 6.9mm, Insert length 200mm, Outer length 190mm, 15mm Bottom length	1	440,00 €

### SHIGEMI TUBE SET for DMSO-D6


Article Number	Description	Quantity	Price
DMS-0025	2,5mm Shigemi TUBE SET for DMSO-D6 for BRUKER Outer OD 2.5mm, Insert OD 1.8mm, Insert length 115mm, Outer length 120mm, 6mm Bottom	1	185,00 €
DMS-003	3mm Shigemi TUBE SET for DMSO-D6 for every System Outer OD 3.0mm, Insert OD 2.4mm, Insert length 190mm, Outer length 180mm, 10mm Bottom length	1	165,00 €
DMS-005B	5mm Shigemi TUBE SET for DMSO-D6 for BRUKER Outer OD 5.0mm, Insert OD 4.1mm, Insert length 190mm, Outer length 180mm, 8mm Bottom length	1	130,00 €
DMS-005J	5mm Shigemi TUBE SET for DMSO-D6 for JOEL Outer OD 5.0mm, Insert OD 4.1mm, Insert length 190mm, Outer length 180mm, 12mm Bottom length	1	155,00 €
DMS-005V	5mm Shigemi TUBE SET for DMSO-D6 for VARIAN Outer OD 5.0mm, Insert OD 4.1mm, Insert length 190mm, Outer length 180mm, 15mm Bottom length	1	155,00 €
DMS-008BJ	8mm Shigemi TUBE SET for DMSO-D6 for BRUKER and JOEL Outer OD 8.0mm, Insert OD 6.9mm, Insert length 200mm, Outer length 190mm, 10mm Bottom length	1	440,00 €
DMS-008V	8mm Shigemi TUBE SET for DMSO-D6 for VARIAN Outer OD 8.0mm, Insert OD 6.9mm, Insert length 200mm, Outer length 190mm, 15mm Bottom length	1	440,00 €
DMS-010BJ	10mm Shigemi TUBE SET for DMSO-D6 for BRUKER and JOEL Outer OD 10.0mm, Insert OD 8.9mm, Insert length 200mm, Outer length 190mm, 10mm Bottom length	1	505,00 €
DMS-010V	10mm Shigemi TUBE SET for DMSO-D6 for VARIAN Outer OD 10.0mm, Insert OD 8.9mm, Insert length 200mm, Outer length 190mm, 15mm Bottom length	1	545,00 €
DMS-005TB	5mm Shigemi Thin Wall TUBE SET ID 4.52mm for DMSO-D6 for BRUKER Outer OD 5mm, Insert OD 4.1mm, Insert length 190mm, Outer length 180mm, 8mm Bottom length	1	220,00 €
DMS-005TJ	5mm Shigemi Thin Wall ID 4.52mm TUBE SET for DMSO-D6 for JOEL Outer OD 5mm, Insert OD 4.1mm, Insert length 190mm, Outer length 180mm, 12mm Bottom length	1	220,00 €
DMS-005TV	5mm Shigemi Thin Wall ID 4.52mm TUBE SET for DMSO-D6 for VARIAN Outer OD 5mm, Insert OD 4.1mm, Insert length 190mm, Outer length 180mm, 15mm Bottom length	1	185,00 €
DMS-008TB	8mm Shigemi Thin Wall TUBE SET ID 7.5mm for D2O for BRUKER Outer OD 8mm, Insert OD 6.9mm, Insert length 200mm, Outer length 190mm, 10mm Bottom length	1	185,00 €
DMS-008TV	8mm Shigemi Thin Wall TUBE SET ID 7.5mm for D2O for VARIAN Outer OD 8mm, Insert OD 6.9mm, Insert length 200mm, Outer length 190mm, 15mm Bottom length	1	465,00 €

### SHIGEMI TUBE SET for Methanol-D4

Article Number	Description	Quantity	Price
MMS-0025	2,5mm Shigemi TUBE SET for Methanol-D4 for BRUKER Outer OD 2.5mm, Insert OD 1.8mm, Insert length 115mm, Outer length 120mm, 6mm Bottom length	1	185,00 €
MMS-003	3mm Shigemi TUBE SET for Methanol-D4 for every System Outer OD 3.0mm, Insert OD 2.4mm, Insert length 190mm, Outer length 180mm, 10mm Bottom length	1	170,00 €
MDMS-005B	5mm Shigemi TUBE SET for Methanol-D4 for BRUKER Outer OD 5.0mm, Insert OD 4.1mm, Insert length 190mm, Outer length 180mm, 8mm Bottom length	1	130,00 €
MMS-005J	5mm Shigemi TUBE SET for Methanol-D4 for JOEL Outer OD 5.0mm, Insert OD 4.1mm, Insert length 190mm, Outer length 180mm, 12mm Bottom length	1	155,00 €
MMS-005V	5mm Shigemi TUBE SET for Methanol-D4 for VARIAN Outer OD 5.0mm, Insert OD 4.1mm, Insert length 190mm, Outer length 180mm, 15mm Bottom length	1	155,00 €
MMS-008BJ	8mm Shigemi TUBE SET for Methanol-D4 for BRUKER and JOEL Outer OD 8.0mm, Insert OD 6.9mm, Insert length 200mm, Outer length 190mm, 10mm Bottom length	1	440,00 €
MMS-008V	8mm Shigemi TUBE SET for Methanol-D4 for VARIAN Outer OD 8.0mm, Insert OD 6.9mm, Insert length 200mm, Outer length 190mm, 15mm Bottom length	1	440,00 €
MMS-010BJ	10mm Shigemi TUBE SET for Methanol-D4 for BRUKER and JOEL Outer OD 10.0mm, Insert OD 8.9mm, Insert length 200mm, Outer length 190mm, 10mm Bottom length	1	505,00 €
MMS-010V	10mm Shigemi TUBE SET for Methanol-D4 for VARIAN Outer OD 10.0mm, Insert OD 8.9mm, Insert length 200mm, Outer length 190mm, 15mm Bottom length	1	545,00 €
MMS-005TB	5mm Shigemi Thin Wall TUBE SET ID 4.52mm for Methanol-D4 for BRUKER Outer OD 5mm, Insert OD 4.1mm, Insert length 190mm, Outer length 180mm, 8mm Bottom length	1	190,00 €
MMS-005TJ	5mm Shigemi Thin Wall ID 4.52mm TUBE SET for Methanol-D4 for JOEL Outer OD 5mm, Insert OD 4.1mm, Insert length 190mm, Outer length 180mm, 12mm Bottom length	1	190,00 €
MMS-005TV	5mm Shigemi Thin Wall ID 4.52mm TUBE SET for Methanol-D4 for VARIAN Outer OD 5mm, Insert OD 4.1mm, Insert length 190mm, Outer length 180mm, 15mm Bottom length	1	215,00 €
MMS-008TB	8mm Shigemi Thin Wall TUBE SET ID 7.5mm for D2O for BRUKER Outer OD 8mm, Insert OD 6.9mm, Insert length 200mm, Outer length 190mm, 10mm Bottom length	1	185,00 €
MMS-008TV	8mm Shigemi Thin Wall TUBE SET ID 7.5mm for D2O for VARIAN Outer OD 8mm, Insert OD 6.9mm, Insert length 200mm, Outer length 190mm, 15mm Bottom length	1	560,00 €

### SHIGEMI Precision Large Volume NMR Tubes

Article Number	Description	Quantity	Price
S-8-001	OD 7.98mm, ID 7.03mm, length 190mm, Straightness 10µm, Wall Thickness control 15µm, 500MHz	1	52,99 €
S-8-002	OD 7.98mm, ID 7.03mm, length 190mm, Straightness 10µm, Wall Thickness control 20µm, 400MHz	1	43,91 €
S-10-001	OD 9.98mm, ID 9.03mm, length 190mm, Straightness 10µm, Wall Thickness control 15µm, 500MHz	1	52,99 €
S-10-002	OD 9.98mm, ID 9.03mm, length 190mm, Straightness 10µm, Wall Thickness control 20µm, 400MHz	1	43,91 €
S-10-003	OD 9.98mm, ID 9.03mm, length 190mm, Straightness 20µm, Wall Thickness control 30µm, 400MHz	1	34,82 €
S-10-004	OD 9.98mm, ID 9.03mm, length 190mm, Straightness 20µm, Wall Thickness control 40µm, 300MHz	1	31,80 €



### SHIGEMI Precision Large Volume NMR Tubes


Article Number	Description	Quantity	Price
PST-001	OD 5mm, Wall Thickness 0.21mm, length Thin Wall 50mm, 600MHz length 180mm, Straightness 10µm, Wall Thickness control 10µm	1	52,99 €
PST-002	OD 5mm, Wall Thickness 0.21mm, length Thin Wall 50mm, 400MHz length 180mm, Straightness 20µm, Wall Thickness control 25µm	1	34,82 €
ST8-001	OD 8mm, Wall Thickness 0.21mm, length Thin Wall 50mm, 600MHz length 190mm, Straightness 10µm, Wall Thickness control 10µm	1	78,73 €
ST8-002	OD 8mm, Wall Thickness 0.21mm, length Thin Wall 50mm, 400MHz length 190mm, Straightness 20µm, Wall Thickness control 20µm	1	62,08 €
ST10-001	OD 10mm, Wall Thickness 0.21mm, length Thin Wall 50mm, 600MHz length 190mm, Straightness 10µm, Wall Thickness control 10µm	1	78,73 €
ST10-002	OD 10mm, Wall Thickness 0.21mm, length Thin Wall 50mm, 400MHz length 190mm, Straightness 20µm, Wall Thickness control 20µm	1	62,08 €






## SHIGEMI Precision Coaxial System

Article Number	Description	Quantity	Price
SP-401-SET	Thin Wall 0,35mm, 183mm length, Upper O.D. innertube 4.1mm, Tube 56ul, Outer Tube 579uL, Stem length 50mm, Stem O.D 1.7mm, Stem I.D. 1.2mm, 1x Inner Tube, 1x Outer Tube, 1x Cap	1	102,96 €
SP-402-SET	Thin Wall 0,35mm, 183mm length, Upper O.D. innertube 4.1mm, Tube 88ul, Outer Tube 535uL Stem length 50mm, Stem O.D 2mm, Stem I.D. 1.5mm, 1x Inner Tube, 1x Outer Tube, 1x Cap	1	102,96 €
SP-403-SET	Thin Wall 0,35mm, 183mm length, Upper O.D. innertube 4.1mm, Tube 141ul, Outer Tube 448uL Stem length 50mm, Stem O.D 2.5mm, Stem I.D. 1.9mm, 1x Inner Tube, 1x Outer Tube, 1x Cap	1	84,79 €
SP-404-SET	Thin Wall 0,35mm, 183mm length, Upper O.D. innertube 4.1mm, Tube 244ul, Outer Tube 341uL Stem length 50mm, Stem O.D 3mm, Stem I.D. 2.5mm, 1x Inner Tube, 1x Outer Tube, 1x Cap	1	84,79 €
SP-405-SET	Thin Wall 0,35mm, 183mm length, Upper O.D. innertube 4.1mm, Tube 399ul, Outer Tube 33uL Stem length 50mm, Stem O.D 4.1mm, Stem I.D. 3.2mm, 1x Inner Tube, 1x Outer Tube, 1x Cap	1	84,79 €
PS-001-7	Thin Wall Outer Tube SP-401/402/403/404/405	1	22,71 €
SPT-401-SET	Ultra Thin Wall 0,21mm, 183mm length, Upper O.D. innertube 4.1mm, Tube 56ul, Outer Tube 771uL Stem length 50mm, Stem O.D 1.7mm, Stem I.D. 1.2mm, 1x Inner Tube, 1x Outer Tube, 1x Cap	1	119,61 €
SPT-402-SET	Ultra Thin Wall 0,21mm, 183mm length, Upper O.D. innertube 4.1mm, Tube 88ul, Outer Tube 728uL Stem length 50mm, Stem O.D 2mm, Stem I.D. 1.5mm, 1x Inner Tube, 1x Outer Tube, 1x Cap	1	119,61 €
SPT-403-SET	Ultra Thin Wall 0,21mm, 183mm length, Upper O.D. innertube 4.1mm, Tube 141ul, Outer Tube 640uL Stem length 50mm, Stem O.D 2.5mm, Stem I.D. 1.9mm, 1x Inner Tube, 1x Outer Tube, 1x Cap	1	102,96 €
SPT-404-SET	Ultra Thin Wall 0,21mm, 183mm length, Upper O.D. innertube 4.1mm, Tube 244ul, Outer Tube 534uL Stem length 50mm, Stem O.D 3mm, Stem I.D. 2.5mm, 1x Inner Tube, 1x Outer Tube, 1x Cap	1	102,96 €
SPT-405-SET	Ultra Thin Wall 0,21mm, 183mm length, Upper O.D. innertube 4.1mm, Tube 399ul, Outer Tube 263uL Stem length 50mm, Stem O.D 4.1mm, Stem I.D. 3.2mm, 1x Inner Tube, 1x Outer Tube, 1x Cap	1	102,96 €
PST-001-7	Thin Wall Outer Tube SP-401/402/403/404/405	1	52,99 €



### SHIGEMI 5mm Micro Bottom Tubes

Article Number	Description	Quantity	Price
SP-501	Brosilicate glass ASTM Type1 Class A, 180mm length, Upper O.D. 5mm Stem length 60mm, Stem O.D 1.7mm, Stem I.D. 1.2mm	1	62,08 €
SP-502	Brosilicate glass ASTM Type1 Class A, 180mm length, Upper O.D. 5m Stem length 60mm, Stem O.D 2mm, Stem I.D. 1.5mm	1	62,08 €
SP-503	Brosilicate glass ASTM Type1 Class A, 180mm length, Upper O.D. 5mm Stem length 60mm, Stem O.D 2.5mm, Stem I.D. 1.9mm	1	43,91 €
SP-504	Brosilicate glass ASTM Type1 Class A, 180mm length, Upper O.D. 5mm Stem length 60mm, Stem O.D 3mm, Stem I.D. 2.5mm	1	43,91 €
SP-505	Brosilicate glass ASTM Type1 Class A, 180mm length, Upper O.D. 5mm Stem length 60mm, Stem O.D 4mm, Stem I.D. 3.2mm	1	43,91 €



### SHIGEMI Quartz NMR Tubes

Article Number	Description	Quantity	Price
SS-002	Shigemi 4.965mm Quartz NMR Tube Fused Silica 190mm length, O.D. 4.965mm, I.D. 4.2mm	1	39,37 €
SS-0010	Shigemi 9.98mm Quartz NMR Tube Fused Silica 190mm length, O.D. 9.98mm, I.D. 9.0mm	1	69,65 €

### SHIGEMI Quartz ESR (EPR) Tubes

Article Number	Description	Quantity	Price
EPS-004J	Shigemi 5mm Quartz ESR (EPR) Tube Fused Silica for BRUKER 270mm length, O.D. 4mm, I.D. 3mm	1	34,82 €
EPS-004B	Shigemi 5mm Quartz ESR (EPR) Tube Fused Silica for JOEL 250mm length, O.D. 4mm, I.D. 3mm	1	34,82 €
EPS-005J	Shigemi 5mm Quartz ESR (EPR) Tube Fused Silica for BRUKER 270mm length, O.D. 5mm, I.D. 4.4mm	1	52,99 €
EPS-005B	Shigemi 5mm Quartz ESR (EPR) Tube Fused Silica 250mm length, O.D. 4mm, I.D. 3mm	1	52,99 €

### SHIGEMI Slot Tube

Article Number	Description	Quantity	Price
SBMS-1	Shigemi 5mm Slot NMR Tube for High concentration electrolyte samples For more information just asks the Rototec-Spintec Team! Do not use Air-lift function!	1	275,00 €
SBMS-5	5x Shigemi 5mm Slot NMR Tube for High concentration electrolyte samples For more information just asks the Rototec-Spintec Team! Do not use Air-lift function!	1	1200,00 €
Positioning-tool-1-Set	For fixing exact position of Slot Tube if no "Bruker Shaped Shigemi Tube Shuttle"	1	2271,15 €

### DURAN® 5mm Ultra Precision Thin Wall NMR Sample Tubes

Article Number	Description	Quantity	Price
PS-001-7	OD 4.965mm, ID 4.22mm, Wall Thickness 0,3725mm, length 180mm Straightness 10µm, Wall Thickness control 5µm, 600MHz	1	33,31 €
PS-001-8	OD 4.965mm, ID 4.22mm, Wall Thickness 0,3725mm, length 200mm Straightness 10µm, Wall Thickness control 5µm, 600MHz	1	33,31 €
PS-002-7	OD 4.965mm, ID 4.22mm, Wall Thickness 0,3725mm, length 180mm Straightness 15µm, Wall Thickness control 10µm, 500MHz	1	28,77 €
PS-002-8	OD 4.965mm, ID 4.22mm, Wall Thickness 0,3725mm, length 200mm Straightness 15µm, Wall Thickness control 10µm, 500MHz	1	30,28 €
PS-003-7	OD 4.965mm, ID 4.22mm, Wall Thickness 0,3725mm, length 180mm Straightness 15µm, Wall Thickness control 20µm, 500MHz	1	25,74 €
PS-003-8	OD 4.965mm, ID 4.22mm, Wall Thickness 0,3725mm, length 200mm Straightness 15µm, Wall Thickness control 20µm, 500MHz	1	27,25 €
PS-004-7	OD 4.965mm, ID 4.22mm, Wall Thickness 0,3725mm, length 180mm Straightness 20µm, Wall Thickness control 25µm, 400MHz	1	16,66 €
PS-004-8	OD 4.965mm, ID 4.22mm, Wall Thickness 0,3725mm, length 200mm Straightness 20µm, Wall Thickness control 25µm, 400MHz	1	18,17 €

### Pyrex® 4mm Ultra Precision Thin Wall NMR Sample Tubes

Article Number	Description	Quantity	Price
PR-001	OD 4mm, ID 3.22mm, length 180mm Straightness 10µm, Wall Thickness control 5µm	1	34,82 €
PR-002	OD 4mm, ID 3.22mm, length 180mm Straightness 10µm, Wall Thickness control 5µm	1	27,25 €

### PYREX® 3mm Ultra Precision Thin Wall NMR Sample Tubes

Article Number	Description	Quantity	Price
PN-001	OD 3mm, ID 2.5mm, Wall Thickness 0,25mm, length 180mm Straightness 15µm, Wall Thickness control 5µm	1	34,82 €
PN-002	OD 3mm, ID 2.5mm, Wall Thickness 0,25mm, length 180mm Straightness 25µm, Wall Thickness control 20µm	1	27,25 €

## Our Partners

Co-Axial Dynamics – Precision RF components for the NMR / MRI community

Daedalus Innovations - Extreme High Pressure Systems for NMR up to and exceeding 3 KBAR

Doty Scientific – Solids, Liquids & MRImaging Probes and accessories

FTS - NMR Sample Cooling & Heating Systems - Plug & Play

Shigemi – Susceptibility Matched NMR Tubes

Tomco Technologies – Linear RF Power Amplifiers for NMR and MRI

TrewMac Vector Impedance analysers ideal for Probe Tuning and other RF applications.

Wilmad LabGlass – Simply the best NMR / EPR Tubes in the world

ZeoChem NMR Deuterated Solvents from Switzerland - Chemistry. Pure. Efficient.

Chamber of Commerce: Stuttgart HRB 762724

Managing Director: Dipl. Ing. David Cross

VAT ID-Nr. DE146616255

E&OE (Errors and Omissions Excluded).

Pricing & Specifications subject to change without notice.

Rototec-Spintec GmbH Terms & Conditions of Sale (German):

[https://www.rototec-spintec.com/docs/AGB\\_GER\\_2017.pdf](https://www.rototec-spintec.com/docs/AGB_GER_2017.pdf)